

eFilm Tools – Quick Reference Guide

Note: The small “R” or “L” on the top of the tool icon indicates that a right or left mouse click or scroll is needed to activate the tool.

To scroll depress and hold the right or left side of the mouse. Scroll by moving the mouse to achieve the desired affect on the image.

Common Tools



Brightness/Contrast Level
Adjusts the brightness and contrast of an image. Right scroll the mouse up to decrease brightness, down to increase brightness. Left to increase contrast, right to increase latitude. Press F2 to return to original brightness/contrast settings



Zoom
Right scroll the mouse up to zoom the image in or out. Scrolling up increases the image zoom. Scrolling down decreases the image zoom.



Panning
Repositions the image in the window. Left scroll the mouse to move the image once it is zoomed. (Use in conjunction with the Zoom tool.)



Reset Magnification
Left click to return the image to the original magnification.



Magnify
Left click on the image to magnify a specific area of the image. Click and drag the mouse over the area of the image you want to magnify. The magnifying window appears and follows the cursor as it magnifies.



Screen Layout
Left click to display the screen layout dialog box. The ‘series of images’ layout determines the format of the panes in the window. Select desired layout.



Toggle Survey / Explode Mode
‘Explodes’ image to fill the screen and returns to the former layout.



Thumbnails
Displays a ‘thumbnail’ image of all images in an open study. The Thumbnail panel also displays a tab for each patient study



Toggle Overlay
Hides or displays the overlay study information and scale bar.



a VCA ANTECH company

Technical Support – 800.268.5354, option 3
For after hours emergencies use the same number, option 3, and select option 1

Image Manipulation Tools



Rotation
Rotates image 90° Counterclockwise



Rotation
Rotates image 90° Clockwise



Flips Image
Flips image left/right or right/left



Flips Image
Flips image top/bottom or bottom/top



Invert
Inverts an image to black on white, or white on black.

Next / Previous Tools



Se Arrows
Selects the ‘previous’ or ‘next’ series in the selected exam.



St Arrows
Selects the ‘previous’ or ‘next’ study from the local exam list.



Measurement Tools



Arrow
To draw an arrow, position cursor at the source (the arrow tail,) right scroll to the destination, the arrow head.) Enter note if desired.



Linear Measurement
Position cursor at starting location, click and right scroll to the ending location. To stretch line: left scroll either end of the line to a new location. To move line: left click anywhere on the line except at the ends and drag it to a new location.



Calibration
Image may be calibrated for magnification factor.



Angle
Left click this tool to measure an angle between two intersecting lines. (To draw intersecting lines, use the Linear Measurement Tool.)



Clear Measurement
Clears ALL arrows, measurements, and angles on current image. To clear a single measurement: select measurement to be cleared, right click and select delete.